

Digital Refractometer for a Water Solution of Hydrogen Peroxide PR-50HO



The concentration of a water solution of hydrogen peroxide used as a disinfecting solution for cartons for soft drinks and milk-based drinks can be measured with only a small quantity of sample and with a simple procedure. The temperature is automatically corrected and not necessary to worry about the temperature.

a	• •	4.
- 51	necific	ations
_	pecine	titions.

Model	PR-50HO
Cat.No.	3478
Scale	Concentration of a water solution of hydrogen peroxide
Range	0.0 to 50.0% (W/W)
Resolution	0.1%
Accuracy	±0.5%
Temperature compensation	5 to 40°C
Power Supply	006P dry battery (9V)
Dimensions & Weight	17×9×4cm, 270g

Options



^{*} Calibration Certificate : Contact an ATAGO representative for details.